



ENCPC -Skin Reliever

Gel Sheath



US Patent
7,709,570

The ALPS Skin Reliever essentially is a gel sheath used to compensate for volume loss. The Skin Reliever is made from ALPS EasyGel, with 3mm of gel at the distal end to help provide distal comfort, tapering down to 1.5mm of gel on the proximal sides to compensate for shrinkage. The Skin Relievers come in three sizes, with the EasyGel having 2000% elasticity, the three sizes will accommodate all patients. See the chart below for sizing.

Skin Reliever ENCP sizing chart

Measurements taken 6cm from the distal end

Measured Size	Part #
10-19 cm	ENCP-16
20 -27cm	ENCP-20
28-37cm	ENCP-28

Donning the Skin Reliever

Although it is not necessary to apply any type of lubricant to the residual limb, you may use either talcum powder or a very small amount of Alps SL108 silicone skin lotion. **Avoid petroleum jelly and other oil based lubricants** as they will damage the gel. **Talcum powder will form a scaling on the surface of the Skin Reliever.** (Using the SL108 is the preferred procedure).

You may don the Skin Reliever by either inverting and rolling it up or donning it like a sock.

Remove any air pockets trapped between the residual limb and the sheath.

Storing the Skin Reliever

Always keep the foam insert in the Skin Reliever when not in use. Laying flat for long periods may cause creases to develop.

Cleaning the Skin Reliever

Daily: Wash the Skin Reliever and the residual limb with warm water and a neutral soap such as ALPS Cleanser PD595 or any other hypoallergenic ph-balanced soap. Use a wash cloth, but do not rub, nor scrub or you will damage the Skin Reliever. Rinse Well. Wipe the Skin Reliever dry with a clean, dry towel.

Weekly: Place a small amount of rubbing alcohol on a soft, clean cloth. Wipe the Skin Reliever with the cloth thoroughly, wipe dry with a towel and apply a light dusting of talcum powder to the outside of the liner. **DO NOT** soak the liner in alcohol.

WARRANTY

30-day satisfaction guaranteed, no questions asked
30-day warranty against manufacturing defects

To obtain an ALPS Return Authorization Number (RA#) call ALPS Customer Service between the hours of 8 am and 5 pm EST and provide an ALPS representative with the following information:

Customer ID # -- Invoice # -- Date of Invoice – Nature of Return

The RA# must be displayed on the exterior of the returned item box or it will be refused at the dock.

Warning: Specific exclusions from this warranty are: cuts, holes, rip due to hits, scrapes, etc.



CEpartner4U
Esdoornlaan 13, 3951DB Maarn, The Netherlands
www.cepartner4u.com

900029 Rev J, CAN# 6653
December 2016
ENCP

ALPS South LLC
2895 42nd Ave. N.
St. Petersburg, FL. 33714
USA
Tel: 1-727-528-8566
Fax: 1-727-528-8862
Info@easyliner.com
www.easyliner.com