


# CUFFIE PROTESICHE

## TABELLA MISURE

### COME SCEGLIERE LA TAGLIA

Misurare la circonferenza a 6 cm dall'apice distale.



CODICE #	DESCRIZIONE	SPESSORE (MM)	TAGLIA										
			16	20	24	26	28	32	38	44			
AKDT AKDTHD AKFR AKFRHD	Extreme Liner	3 o 6											
APDT-TP APFR-TP	Anterior Posterior Tapered Liner	6/3		20	24	26	28	32					
ECDTHD ECFRHD	ECO Liner	3 o 6		20	24	26	28	32	38	44			
ELDT ELFR	EasyLiner	3, 6, 6/3 o 9*	16	20	24*	26	28*	32	38	44			
GPDT GPFR	General Purpose	3 o 6	16	20	24	26	28	32	38	44			
SPDT SPDTHD SPFR SPFRHD	Superior Performance	3 o 6	16	20	24	26	28	32	38	44			
SADT SAFR SBFR	Smart Seal SADT	3 o 6		20	24	26	28	32	38				
	Smart Seal SAFR	3 o 6				26	28	32	38				
	Smart Seal SBFR	3 o 6		20	24	26	28						
EZF	EasyFlex	3 o 6	16	20	24	26	28	32	38	44			
VSDT	Liberty Liner	6		20	24	26	28	32	38				
PSI	Proximal Seal Interface	6 o 6/3	16	20	24	26	28	32	38	44			
CIRCONFERENZA (CM)			16-19	20-23	24-25	26-27	28-31	32-37	38-44	44-53			

\*9mm disponibile solo nelle taglie 24 e 28

CODICE #	DESCRIZIONE	SPESSORE (MM)	TAGLIA											
			20	22	24	26	28	30	32	35	38	41	44	
LSC LSL	Silicone Pro	3												
CIRCONFERENZA (CM)			20-22	22-24	24-26	26-28	28-30	30-32	32-35	35-38	38-41	41-44	44-47	

CODICE #	DESCRIZIONE	SPESSORE (MM)	TAGLIA									
			10	16	20	24	26	28	32	38		
ELPX	EasyLiner SuperStretch	6										
CIRCONFERENZA (CM)			10-15	16-19	20-23	24-25	26-27	28-31	32-37	38-43		

**AALPS**

Tel +39 051 86 06 42 | Fax +39 051 86 59 668

www.easyliner.eu | info@alpsitaly.com

V.le della Mercanzia, 85 | 40050 Argelato (BO) - Italy