


ULTRASEAL SILICONE LINER

A step above.

USHSL | USHSC | USLSL | USLSC | USUSL | USUSC



DESCRIPTION

Experience the innovative UltraSeal silicone liner, carefully crafted to meet the needs of both transtibial and transfemoral amputees at all activity levels.

With an integrated band available at three different heights, users will enjoy an airtight seal upon socket entry. Combined with OptiSil, it provides exceptional comfort and optimized suspension for maximum convenience in daily use.

Additionally, in the locking version, the UltraSeal offers a dual suspension system that seamlessly harnesses the advantages of both locking suspension and Sealing Technology, delivering unmatched functionality.

PART #	SUSPENSION
USHSL/USLSL/USUSL	LOCKING
USHSC/USLSC/USUSC	CUSHION

Recommended for all activity levels.



WHY ULTRA?

1

LESS CONSTRICTIVE.

Outer premium fabric combines strength, durability, and differentiated stretch for comfortable elasticity.

2

LOCKING SECURITY.

The locking feature ensures users are provided with a suspension that is secure, firm, and comfortable.


3

FRIENDLY SILICONE SURFACE.


OptiSil offers a thin, lightweight, and smooth inner surface with 1/5 the friction of regular silicone, ensuring superior skin conformance.

CROSS SECTION DIAGRAM


1 THICKNESS | 3 RING HEIGHT OPTIONS | CUSHION & LOCKING



HIGH ring begins 18 cm
from the center of the distal end.




LOW ring begins 13 cm
from the center of the distal end.









ULTRA LOW ring begins 8 cm
from the center of the distal end.

HOW TO SIZE

Measurements taken 6cm
from the distal end



SIZING CHARTS

MEASURED SIZE (CM)	 HIGH RING LOCKING 2MM UNIFORM	 HIGH RING CUSHION 2MM UNIFORM	 LOW RING LOCKING 2MM UNIFORM	 LOW RING CUSHION 2MM UNIFORM	 ULTRA LOW RING LOCKING 2MM UNIFORM	 ULTRA LOW RING CUSHION 2MM UNIFORM
20-21 cm	USHSL20	USHSC20	USLSL20	USLSC20	USUSL20	USUSC20
22-23 cm	USHSL22	USHSC22	USLSL22	USLSC22	USUSL22	USUSC22
24-25 cm	USHSL24	USHSC24	USLSL24	USLSC24	USUSL24	USUSC24
26-27 cm	USHSL26	USHSC26	USLSL26	USLSC26	USUSL26	USUSC26
28-29 cm	USHSL28	USHSC28	USLSL28	USLSC28	USUSL28	USUSC28
30-31 cm	USHSL30	USHSC30	USLSL30	USLSC30	USUSL30	USUSC30
32-34 cm	USHSL32	USHSC32	USLSL32	USLSC32	USUSL32	USUSC32
35-37 cm	USHSL35	USHSC35	USLSL35	USLSC35	USUSL35	USUSC35
38-40 cm	USHSL38	USHSC38	USLSL38	USLSC38	USUSL38	USUSC38
41-43 cm	USHSL41	USHSC41	USLSL41	USLSC41	USUSL41	USUSC41
44-47 cm	USHSL44	USHSC44	USLSL44	USLSC44	USUSL44	USUSC44