

# ANTERIOR POSTERIOR TAPERED

## APDT-TP | APFR-TP


Gel Liner —  
With or Without Distal Attachment



### DESCRIPTION

The ALPS™ AP Tapered Liner is gradually tapered from the anterior to the posterior to provide superior comfort. This liner is available in a pin and lock system or suction suspension. The AP Tapered Liner features our black high performance fabric with gel to assist in reducing bunching in the popliteal region.

### TOP FEATURES



▶ Constructed using high performance black fabric with gel to reduce bunching in the popliteal region

▶ Tapered - 6mm Anterior to 3mm Posterior

MAKING LIVES BETTER




#makinglivesbetter




# CROSS SECTION DIAGRAM

6.0 MM TO 3.0 MM TAPERED



## HOW TO SIZE

Measurements taken 6cm from the distal end



## SIZING CHART

MEASURED SIZE (CM)	LOCKING 6/3 MM TAPERED	CUSHION 6/3 MM TAPERED
20-23 cm	APDT <b>20</b> - TP	APFR <b>20</b> - TP
24-25 cm	APDT <b>24</b> - TP	APFR <b>24</b> - TP
26-27 cm	APDT <b>26</b> - TP	APFR <b>26</b> - TP
28-31 cm	APDT <b>28</b> - TP	APFR <b>28</b> - TP
32-37 cm	APDT <b>32</b> - TP	APFR <b>32</b> - TP



Tel +420 377 223 127 | Fax +420 377 223 010

[www.easyliner.eu](http://www.easyliner.eu) | [info@easyliner.eu](mailto:info@easyliner.eu)

Božkovské náměstí 17/21, 32600 Plzeň - Czech Republic

MAKING  
LIVES  
BETTER